

PH #1066 Sterilization

Application: **STERILIZATION**

The customer requested a Heat Exchanger to be installed on a cold PW line at 20C. The loop is heated once a week with steam to 80C for **Sterilization** and is cooled back to 20C using water.

Solution

- Sanitary Shell & Tube Heat Exchanger
- Model 02125, 95 series
- All components 316L Stainless Steel

Features & Benefits

- Heat Transfer Area: 33.79ft²
- Double Tubesheet
- Shell Diameter: 4"
- Tube Length: 31"
- Tube Count: 421



How to contact Exergy

Please call 1-516-832-9300
or visit www.exergyllc.com
to locate your local rep.